

兔抗 AGRP 多克隆抗体

中文名称：兔抗 AGRP 多克隆抗体

英文名称： Anti-AGRP rabbit polyclonal antibody

别名： ART; AGRT; ASIP2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： AGRP

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an antagonist of the melanocortin-3 and melanocortin-4 receptor. It appears to regulate hypothalamic control of feeding behavior via melanocortin receptor and/or intracellular calcium regulation, and thus plays a role in weight homeostasis. Mutations in this gene have been associated with late on-set obesity.
Applications:	ELISA, IHC
Name of antibody:	AGRP
Immunogen:	Synthetic peptide of human AGRP
Full name:	agouti related neuropeptide

Synonyms:	ART; AGRT; ASIP2
SwissProt:	O00253
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human gastric cancer
IHC Recommend dilution:	25-100

