

兔抗 AGTRAP 多克隆抗体

中文名称：兔抗 AGTRAP 多克隆抗体

英文名称：Anti-AGTRAP rabbit polyclonal antibody

别名：ATRAP

相关类别：一抗

储存：冷冻（-20℃）

抗原：AGTRAP

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	AGTRAP
Immunogen:	Fusion protein of human AGTRAP
Full name:	angiotensin II receptor-associated protein
Synonyms:	ATRAP

SwissProt:	Q6RW13
IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	50-200

