

## 兔抗 AHR(Phospho-Ser36) 多克隆抗体

中文名称：兔抗 AHR(Phospho-Ser36) 多克隆抗体

英文名称： Anti-AHR(Phospho-Ser36) rabbit polyclonal antibody

别名： bHLHe76

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： AHR(Phospho-Ser36)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a ligand-activated transcription factor involved in the regulation of biological responses to planar aromatic hydrocarbons. This receptor has been shown to regulate xenobiotic-metabolizing enzymes such as cytochrome P450. Its ligands included a variety of aromatic hydrocarbons.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	AHR(Phospho-Ser36)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Serine 36 (N-P-S(p)-K-R) derived from Human AhR.
<b>Full name:</b>	aryl hydrocarbon receptor

<b>Synonyms :</b>	bHLHe76
<b>SwissProt:</b>	P35869
<b>WB Predicted band size:</b>	96 kDa
<b>WB Positive control:</b>	HepG2 cells lysates
<b>WB Recommended dilution:</b>	500-1000

