

兔抗 AHR 多克隆抗体

中文名称：兔抗 AHR 多克隆抗体

英文名称：Anti-AHR rabbit polyclonal antibody

别名：bHLHe76

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AHR

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a ligand-activated transcription factor involved in the regulation of biological responses to planar aromatic hydrocarbons. This receptor has been shown to regulate xenobiotic-metabolizing enzymes such as cytochrome P450. Its ligands included a variety of aromatic hydrocarbons.
Applications:	ELISA, IHC
Name of antibody:	AHR
Immunogen:	Synthetic peptide of human AHR
Full name:	aryl hydrocarbon receptor
Synonyms :	bHLHe76

SwissProt:	P35869
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human cervical cancer
IHC Recommend dilution:	50-200

