

兔抗 CHD1L 多克隆抗体

中文名称：兔抗 CHD1L 多克隆抗体

英文名称：Anti-CHD1L rabbit polyclonal antibody

别名：chromodomain helicase DNA binding protein 1 like; ALC1; CHDL

相关类别：一抗

抗原：CHD1L

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a DNA helicase protein involved in DNA repair. The protein converts ATP to add poly(ADP-ribose) as it regulates chromatin relaxation following DNA damage. Several alternatively spliced transcripts variants have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	CHD1L
Immunogen:	Synthetic peptide of human CHD1L
Full name:	chromodomain helicase DNA binding protein 1 like
Synonyms:	ALC1; CHDL
SwissProt:	Q86WJ1

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100

