

兔抗 AKAP6 多克隆抗体

中文名称：兔抗 AKAP6 多克隆抗体

英文名称： Anti-AKAP6 rabbit polyclonal antibody

别名： ADAP6; PRKA6; mAKAP; ADAP100; AKAP100

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： AKAP6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The A-kinase anchor proteins (AKAPs) are a group of structurally diverse proteins, which have the common function of binding to the regulatory subunit of protein kinase A (PKA) and confining the holoenzyme to discrete locations within the cell. This gene encodes a member of the AKAP family. The encoded protein is highly expressed in various brain regions and cardiac and skeletal muscle. It is specifically localized to the sarcoplasmic reticulum and nuclear membrane, and is involved in anchoring PKA to the nuclear membrane or sarcoplasmic reticulum.

Applications:	ELISA, IHC
Name of antibody:	AKAP6
Immunogen:	Synthetic peptide of human AKAP6
Full name:	A kinase (PRKA) anchor protein 6
Synonyms :	ADAP6; PRKA6; mAKAP; ADAP100; AKAP100
SwissProt:	Q13023
ELISA Recommended dilution:	3000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

