

兔抗 CHDH 多克隆抗体

中文名称：兔抗 CHDH 多克隆抗体

英文名称： Anti-CHDH rabbit polyclonal antibody

相关类别： 一抗

抗 原： CHDH

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a choline dehydrogenase that localizes to the mitochondrion. Variations in this gene can affect susceptibility to choline deficiency. A few transcript variants have been found for this gene, but the full-length nature of only one has been characterized to date.
Applications:	ELISA, WB, IHC
Name of antibody:	CHDH
Immunogen:	Fusion protein of human CHDH
Full name:	choline dehydrogenase
SwissProt:	Q8NE62
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human breast cancer

IHC Recommend dilution:	25-100
WB Predicted band size:	65 kDa
WB Positive control:	231 cells, mouse kidney and liver tissue
WB Recommended dilution:	500-2000



