

兔抗 AKR1B10 多克隆抗体

中文名称：兔抗 AKR1B10 多克隆抗体

英文名称：Anti-AKR1B10 rabbit polyclonal antibody

别名：HIS; HSI; ARL1; ARL-1; ALDRLn; AKR1B11; AKR1B12

相关类别：一抗

储存：冷冻 (-20℃)

抗原：AKR1B10

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the aldo/keto reductase superfamily, which consists of more than 40 known enzymes and proteins. This member can efficiently reduce aliphatic and aromatic aldehydes, and it is less active on hexoses. It is highly expressed in adrenal gland, small intestine, and colon, and may play an important role in liver carcinogenesis.
Applications:	ELISA, WB, IHC
Name of antibody:	AKR1B10
Immunogen:	Fusion protein of human AKR1B10
Full name:	aldo-keto reductase family 1, member B10 (aldose

	reductase)
Synonyms :	HIS; HSI; ARL1; ARL-1; ALDRLn; AKR1B11; AKR1B12
SwissProt:	O60218
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	36 kDa
WB Positive control:	Raji and Jurkat cells
WB Recommended dilution:	200-1000



