

兔抗 CHI3L2 多克隆抗体

中文名称：兔抗 CHI3L2 多克隆抗体

英文名称： Anti-CHI3L2 rabbit polyclonal antibody

别名： CHIL2; YKL39; YKL-39

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CHI3L2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is similar to bacterial chitinases but lacks chitinase activity. The encoded protein is secreted and is involved in cartilage biogenesis. Several transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	CHI3L2
Immunogen:	Fusion protein corresponding to a region derived from 191-390 amino acids of human CHI3L2
Full name:	chitinase 3 like 2
Synonyms:	CHIL2; YKL39; YKL-39

SwissProt:	Q15782
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human liver cancer;Human tonsil
IHC Recommended dilution:	25-100



