

兔抗 CHML 多克隆抗体

中文名称: 兔抗 CHML 多克隆抗体

英文名称: Anti-CHML rabbit polyclonal antibody

别名: REP2

抗原: CHML

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	he product of the CHML gene supports geranylgeranylation of most Rab proteins and may substitute for REP-1 in tissues other than retina. CHML is localized close to the gene for Usher syndrome type II.
Applications:	WB, IHC
Name of antibody:	CHML
Immunogen:	Synthesized peptide derived from internal of human CHML.
Full name:	choroideremia-like (Rab escort protein 2)
Synonyms :	REP2
SwissProt:	P26374

IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	74 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000



