

## 兔抗 CHMP4B 多克隆抗体

中文名称：兔抗 CHMP4B 多克隆抗体

英文名称： Anti-CHMP4B rabbit polyclonal antibody

别名： SNF7; CTPP3; Shax1; CHMP4A; SNF7-2; VPS32B; CTRCT31; Vps32-2; C20orf178; dJ553F4.4

相关类别： 一抗

抗原： CHMP4B

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the chromatin-modifying protein/charged multivesicular body protein (CHMP) protein family. The protein is part of the endosomal sorting complex required for transport (ESCRT) complex III (ESCRT-III), which functions in the sorting of endocytosed cell-surface receptors into multivesicular endosomes. The ESCRT machinery also functions in the final abscission stage of cytokinesis and in the budding of enveloped viruses such as HIV-1.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	CHMP4B
<b>Immunogen:</b>	Fusion protein of human CHMP4B
<b>Full name:</b>	charged multivesicular body protein 4B
<b>Synonyms :</b>	SNF7; CTPP3; Shax1; CHMP4A; SNF7-2; VPS32B; CTRCT3 1; Vps32-2; C20orf178; dJ553F4.4
<b>SwissProt:</b>	Q9H444
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human liver cancer
<b>IHC Recommend dilution:</b>	50-200

