

## 兔抗 ALDH1B1 多克隆抗体

中文名称：兔抗 ALDH1B1 多克隆抗体

英文名称： Anti-ALDH1B1 rabbit polyclonal antibody

别名： ALDH5; ALDHX

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： ALDH1B1

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This protein belongs to the aldehyde dehydrogenases family of proteins. Aldehyde dehydrogenase is the second enzyme of the major oxidative pathway of alcohol metabolism. This gene does not contain introns in the coding sequence. The variation of this locus may affect the development of alcohol-related problems.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	ALDH1B1
<b>Immunogen:</b>	Synthetic peptide of human ALDH1B1
<b>Full name:</b>	aldehyde dehydrogenase 1 family, member B1

<b>Synonyms :</b>	ALDH5; ALDHX
<b>SwissProt:</b>	P30837
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	57 kDa
<b>WB Positive control:</b>	Human liver cancer tissue
<b>WB Recommended dilution:</b>	500-2000

