

兔抗 CHPF 多克隆抗体

中文名称：兔抗 CHPF 多克隆抗体

英文名称： Anti-CHPF rabbit polyclonal antibody

别名： CSS2; CHSY2

抗原： CHPF

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Has both beta-1,3-glucuronic acid and beta-1,4-N-acetylgalactosamine transferase activity. Transfers glucuronic acid (GlcUA) from UDP-GlcUA and N-acetylgalactosamine (GalNAc) from UDP-GalNAc to the non-reducing end of the elongating chondroitin polymer. Isoform 2 may facilitate PARK2 transport into the mitochondria. In collaboration with PARK2, isoform 2 may enhance cell viability and protect cells from oxidative stress.
Applications:	WB, IHC, IF
Name of antibody:	CHPF
Immunogen:	Synthesized peptide derived from internal of human CHSS

	2.
Full name:	chondroitin polymerizing factor
Synonyms :	CSS2; CHSY2
SwissProt:	Q8IZ52
IHC positive control:	Human cervix carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	85 kDa
WB Positive control:	LOVO cells and K562 cells lysates
WB Recommended dilution:	500-3000
IF positive control:	MCF-7 cells
IF Recommend dilution:	100-500



