

兔抗 CHRAC1 多克隆抗体

中文名称：兔抗 CHRAC1 多克隆抗体

英文名称：Anti-CHRAC1 rabbit polyclonal antibody

别名：chromatin accessibility complex 1; YCL1; CHARC1; CHARC15; CHRAC-1; CHRAC15;
CHRAC-15

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：CHRAC1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	CHRAC1 is a histone-fold protein that interacts with other histone-fold proteins to bind DNA in a sequence-independent manner. These histone-fold protein dimers combine within larger enzymatic complexes for DNA transcription, replication, and packaging.
Applications:	ELISA, IHC
Name of antibody:	CHRAC1
Immunogen:	Fusion protein of human CHRAC1
Full name:	chromatin accessibility complex 1
Synonyms:	YCL1; CHARC1; CHARC15; CHRAC-1; CHRAC15; CHRAC-15

SwissProt:	Q9NRG0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100

