

兔抗 CHRDL1 多克隆抗体

中文名称：兔抗 CHRDL1 多克隆抗体

英文名称： Anti-CHRDL1 rabbit polyclonal antibody

别名： CHL; MGC1; MGCN; VOPT; NRLN1; dA141H5.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CHRDL1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an antagonist of bone morphogenetic protein 4. The encoded protein may play a role in topographic retinotectal projection and in the regulation of retinal angiogenesis in response to hypoxia. Alternatively spliced transcript variants encoding different isoforms have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	CHRDL1
Immunogen:	Synthetic peptide of humanCHRDL1
Full name:	chordin-like 1
Synonyms :	CHL; MGC1; MGCN; VOPT; NRLN1; dA141H5.1

SwissProt:	Q9BU40
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	51 kDa
WB Positive control:	293T cells and Human normal kidney tissue
WB Recommended dilution:	200-1000



