

## 兔抗 CHRM1 多克隆抗体

中文名称：兔抗 CHRM1 多克隆抗体

英文名称： Anti-CHRM1 rabbit polyclonal antibody

别名： M1; HM1; M1R

抗原： CHRM1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

The muscarinic cholinergic receptors belong to a larger family of G protein-coupled receptors. The functional diversity of these receptors is defined by the binding of acetylcholine and includes cellular responses such as adenylate cyclase inhibition, phosphoinositide degeneration, and potassium channel mediation. Muscarinic receptors influence many effects of acetylcholine in the central and peripheral nervous system. The muscarinic cholinergic receptor 1 is involved in mediation of vagally-induced bronchoconstriction and in the acid secretion of the gastrointestinal tract. The gene encoding this receptor is localized to 11q13

<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CHRM1
<b>Immunogen:</b>	Fusion protein of human CHRM1
<b>Full name:</b>	cholinergic receptor, muscarinic 1
<b>Synonyms :</b>	M1; HM1; M1R
<b>SwissProt:</b>	P11229
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	51 kDa
<b>WB Positive control:</b>	Mouse brain tissue
<b>WB Recommended dilution:</b>	200-1000

