

兔抗 ALG6 多克隆抗体

中文名称：兔抗 ALG6 多克隆抗体

英文名称：Anti-ALG6 rabbit polyclonal antibody

别名：CDG1C

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ALG6

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ALG6/ALG8 glucosyltransferase family. The encoded protein catalyzes the addition of the first glucose residue to the growing lipid-linked oligosaccharide precursor of N-linked glycosylation. Mutations in this gene are associated with congenital disorders of glycosylation type Ic.
Applications:	ELISA, WB
Name of antibody:	ALG6
Immunogen:	Synthetic peptide of human ALG6
Full name:	ALG6, alpha-1,3-glucosyltransferase
Synonyms :	CDG1C

SwissProt:	Q9Y672
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	58 kDa
WB Positive control:	Human lymphoma tissue
WB Recommended dilution:	200-1000

