

兔抗 ALG9 多克隆抗体

中文名称：兔抗 ALG9 多克隆抗体

英文名称：Anti-ALG9 rabbit polyclonal antibody

别名：CDG1L; DIBD1; LOH11CR1J

相关类别：一抗

抗原：ALG9

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an alpha-1,2-mannosyltransferase enzyme that functions in lipid-linked oligosaccharide assembly. Mutations in this gene result in congenital disorder of glycosylation type II. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	ALG9
Immunogen:	Fusion protein of human ALG9
Full name:	ALG9, alpha-1,2-mannosyltransferase
Synonyms :	CDG1L; DIBD1; LOH11CR1J
SwissProt:	Q9H6U8

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	69 kDa
WB Positive control:	Human fetal liver tissue
WB Recommended dilution:	200-1000



