

兔抗 CHRNA5 多克隆抗体

中文名称：兔抗 CHRNA5 多克隆抗体

英文名称：Anti-CHRNA5 rabbit polyclonal antibody

别名：LNCR2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CHRNA5

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene is a nicotinic acetylcholine receptor subunit and a member of a superfamily of ligand-gated ion channels that mediate fast signal transmission at synapses. These receptors are thought to be heteropentamers composed of separate but similar subunits. Defects in this gene have been linked to susceptibility to lung cancer type 2 (LNCR2).
Applications:	ELISA, WB, IHC
Name of antibody:	CHRNA5
Immunogen:	Synthetic peptide of human CHRNA5
Full name:	cholinergic receptor, nicotinic, alpha 5 (neuronal)

Synonyms :	LNCR2
SwissProt:	P30532
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	53 kDa
WB Positive control:	Jurkat cell
WB Recommended dilution:	200-1000



