

兔抗 ALK(Phospho-Tyr1096) 多克隆抗体

中文名称：兔抗 ALK(Phospho-Tyr1096) 多克隆抗体

英文名称： Anti-ALK(Phospho-Tyr1096) rabbit polyclonal antibody

别名： CD246; NBLST3

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： ALK(Phospho-Tyr1096)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a receptor tyrosine kinase, which belongs to the insulin receptor superfamily. This protein comprises an extracellular domain, an hydrophobic stretch corresponding to a single pass transmembrane region, and an intracellular kinase domain. It plays an important role in the development of the brain and exerts its effects on specific neurons in the nervous system. This gene has been found to be rearranged, mutated, or amplified in a series of tumours including anaplastic large cell lymphomas, neuroblastoma, and non-small cell lung cancer.
Applications:	IHC

Name of antibody:	ALK(Phospho-Tyr1096)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 1096 (P-N-Y(p)-C-F) derived from Human ALK.
Full name:	anaplastic lymphoma receptor tyrosine kinase
Synonyms :	CD246; NBLST3
SwissProt:	Q9UM73
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

