

兔抗 ALKBH1 多克隆抗体

中文名称：兔抗 ALKBH1 多克隆抗体

英文名称： Anti-ALKBH1 rabbit polyclonal antibody

别名： ABH; ABH1; alkB; hABH; ALKBH

抗原： ALKBH1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a homolog to the E. coli alkB gene product. The E. coli alkB protein is part of the adaptive response mechanism of DNA alkylation damage repair. It is involved in damage reversal by oxidative demethylation of 1-methyladenine and 3-methylcytosine.
Applications:	ELISA, WB
Name of antibody:	ALKBH1
Immunogen:	Fusion protein of human ALKBH1
Full name:	alkB homolog 1, histone H2A dioxygenase
Synonyms:	ABH; ABH1; alkB; hABH; ALKBH

SwissProt:	Q13686
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	44 kDa
WB Positive control:	Mouse pancreas tissue
WB Recommended dilution:	500-2000

