

兔抗 CHRNA9 多克隆抗体

中文名称: 兔抗 CHRNA9 多克隆抗体

英文名称: Anti-CHRNA9 rabbit polyclonal antibody

别 名: NACHRA9; HSA243342

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: CHRNA9

反应种属: Human, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene is a member of the ligand-gated ionic channe I family and nicotinic acetylcholine receptor gene superfamily. It encodes a plasma membrane protein that forms homo- or hetero-oligomeric divalent cation channels. This protein is involved in cochlea hair cell development and is also expressed in the outer hair cells (OHCs) of the adult cochlea.
Applications:	ELISA, WB, IHC
Name of antibody:	CHRNA9
Immunogen:	Synthetic peptide of human CHRNA9
Full name:	cholinergic receptor nicotinic alpha 9 subunit



Synonyms:	NACHRA9; HSA243342
SwissProt:	Q9UGM1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	40-200
WB Predicted band size:	55 kDa
WB Positive control:	Human cerebrum tissue and Human cerebella tissue lysa tes
WB Recommended dilution:	200-1000





