

兔抗 CHRNA10 多克隆抗体

中文名称：兔抗 CHRNA10 多克隆抗体

英文名称： Anti-CHRNA10 rabbit polyclonal antibody

别名： no

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CHRNA10

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Ionotropic receptor with a probable role in the modulation of auditory stimuli. Agonist binding may induce an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane. The channel is permeable to a range of divalent cations including calcium, the influx of which may activate a potassium current which hyperpolarizes the cell membrane. In the ear, this may lead to a reduction in basilar membrane motion, altering the activity of auditory nerve fibers and reducing the range of dynamic hearing. This may protect against acoustic trauma.

	ma.
Applications:	ELISA, IHC
Name of antibody:	CHRNA10
Immunogen:	Synthetic peptide of human CHRNA10
Full name:	cholinergic receptor, nicotinic, alpha 10 (neuronal)
Synonyms :	no
SwissProt:	Q9GZZ6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human lung cancer
IHC Recommend dilution:	25-100

