

兔抗 ALKBH2 多克隆抗体

中文名称：兔抗 ALKBH2 多克隆抗体

英文名称： Anti-ALKBH2 rabbit polyclonal antibody

别名： ABH2

储存： 冷冻（-20℃）

抗原： ALKBH2

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The Escherichia coli AlkB protein protects against the cytotoxicity of methylating agents by repair of the specific DNA lesions generated in single-stranded DNA. ALKBH2 and ALKBH3 (MIM 610603) are E. coli AlkB homologs that catalyze the removal of 1-methyladenine and 3-methylcytosine (Duncan et al., 2002 [PubMed 12486230]).
Applications:	ELISA, WB
Name of antibody:	ALKBH2
Immunogen:	Fusion protein of human ALKBH2
Full name:	alkB homolog 2, alpha-ketoglutarate dependent dioxy

	genase
Synonyms:	ABH2
SwissProt:	Q6NS38
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	29 kDa
WB Positive control:	Human fetal liver and seminoma tissue, RAW264.7 cells
WB Recommended dilution:	500-2000

