

兔抗 CHRNB1 多克隆抗体

中文名称：兔抗 CHRNB1 多克隆抗体

英文名称： Anti-CHRNB1 rabbit polyclonal antibody

别名： ACHRB; CHRNB; CMS1D; CMS2A; SCCMS

抗原： CHRNB1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The muscle acetylcholine receptor is composed of five subunits: two alpha subunits and one beta, one gamma, and one delta subunit. This gene encodes the beta subunit of the acetylcholine receptor. The acetylcholine receptor changes conformation upon acetylcholine binding leading to the opening of an ion-conducting channel across the plasma membrane. Mutations in this gene are associated with slow-channel congenital myasthenic syndrome.
Applications:	WB
Name of antibody:	CHRNB1
Immunogen:	Synthesized peptide derived from N-terminal of human CHR

	NB1.
Full name:	cholinergic receptor, nicotinic, beta 1 (muscle)
Synonyms :	ACHRB; CHRNB; CMS1D; CMS2A; SCCMS
SwissProt:	P11230
WB Predicted band size:	57 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000

