

## 兔抗 ALPG 多克隆抗体

中文名称：兔抗 ALPG 多克隆抗体

英文名称： Anti-ALPG rabbit polyclonal antibody

别名： alkaline phosphatase, germ cell; GCAP; ALPPL; ALPP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ALPG

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The product of this gene is a membrane bound glycosylated enzyme, localized to testis, thymus and certain germ cell tumors, that is closely related to both the placental and intestinal forms of alkaline phosphatase. [provided by RefSeq, Jul 2008]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ALPG
<b>Immunogen:</b>	Synthetic peptide of human ALPG
<b>Full name:</b>	alkaline phosphatase, germ cell

<b>Synonyms:</b>	GCAP; ALPPL; ALPPL2
<b>SwissProt:</b>	P10696
<b>IHC positive control:</b>	Human esophagus cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

