

## 兔抗 ALPK1 多克隆抗体

中文名称: 兔抗 ALPK1 多克隆抗体

英文名称: Anti-ALPK1 rabbit polyclonal antibody

别名: LAK; 8430410J10Rik

相关类别: 一抗

抗原: ALPK1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an alpha kinase. Mice which were homozygous for disrupted copies of this gene exhibited coordination defects (PMID: 21208416). Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ALPK1
<b>Immunogen:</b>	Fusion protein of human ALPK1
<b>Full name:</b>	alpha kinase 1
<b>Synonyms:</b>	LAK; 8430410J10Rik
<b>SwissProt:</b>	Q96QP1

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human brain
<b>IHC Recommend dilution:</b>	50-300



