

兔抗 CHST2 多克隆抗体

中文名称：兔抗 CHST2 多克隆抗体

英文名称：Anti-CHST2 rabbit polyclonal antibody

别名：carbohydrate sulfotransferase 2; C6ST; GST2; GST-2; Gn6ST-1; HEL-S-75; glcNAc6ST-1

抗原：CHST2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This locus encodes a sulfotransferase protein. The encoded enzyme catalyzes the sulfation of a nonreducing N-acetylglucosamine residue, and may play a role in biosynthesis of 6-sulfosialyl Lewis X antigen.
Applications:	ELISA, IHC
Name of antibody:	CHST2
Immunogen:	Synthetic peptide of human CHST2
Full name:	carbohydrate sulfotransferase 2
Synonyms:	C6ST; GST2; GST-2; Gn6ST-1; HEL-S-75; glcNAc6ST-1
SwissProt:	Q9Y4C5

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-50

