

兔抗 CHST3 多克隆抗体

中文名称：兔抗 CHST3 多克隆抗体

英文名称：Anti-CHST3 rabbit polyclonal antibody

别名：carbohydrate sulfotransferase 3; HSD; C6ST; C6ST1

相关类别：一抗

抗原：CHST3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an enzyme which catalyzes the sulfation of chondroitin, a proteoglycan found in the extracellular matrix and most cells which is involved in cell migration and differentiation. Mutations in this gene are associated with spondylepiphyseal dysplasia and humerospinal dysostosis.
Applications:	ELISA, WB, IHC
Name of antibody:	CHST3
Immunogen:	Synthetic peptide of human CHST3
Full name:	carbohydrate sulfotransferase 3
Synonyms:	HSD; C6ST; C6ST1

SwissProt:	Q7LGC8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and Human gastric cancer
IHC Recommend dilution:	50-100
WB Predicted band size:	55 kDa
WB Positive control:	Mouse heart tissue lysates
WB Recommended dilution:	500-2000



