

## 兔抗 CHST10 多克隆抗体

中文名称：兔抗 CHST10 多克隆抗体

英文名称：Anti-CHST10 rabbit polyclonal antibody

别名：carbohydrate sulfotransferase 10; HNK1ST; HNK-1ST

相关类别：一抗

储存：冷冻（-20℃）

抗原：CHST10

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This protein encoded by this gene transfers sulfate to the C-3 hydroxyl of terminal glucuronic acid of protein- and lipid-linked oligosaccharides. This protein was first identified as a sulfotransferase that acts on the human natural killer-1 (HNK-1) glycan; HNK-1 is a carbohydrate involved in neurodevelopment and synaptic plasticity.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CHST10
<b>Immunogen:</b>	Fusion protein of human CHST10
<b>Full name:</b>	carbohydrate sulfotransferase 10
<b>Synonyms:</b>	HNK1ST; HNK-1ST

<b>SwissProt:</b>	O43529
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	PC-3 cell lysate
<b>WB Recommended dilution:</b>	500-2000





