

兔抗 AMER1 多克隆抗体

中文名称：兔抗 AMER1 多克隆抗体

英文名称：Anti-AMER1 rabbit polyclonal antibody

别名：WTX; OSCS; FAM123B; RP11-403E24.2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AMER1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene upregulates transcriptional activation by the Wilms tumor protein and interacts with many other proteins, including CTNNB1, APC, AXIN1, and AXIN2. Defects in this gene are a cause of osteopathia striata with cranial sclerosis (OSCS).
Applications:	ELISA, IHC
Name of antibody:	AMER1
Immunogen:	Synthetic peptide of human AMER1
Full name:	APC membrane recruitment protein 1
Synonyms :	WTX; OSCS; FAM123B; RP11-403E24.2

SwissProt:	Q5JTC6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	50-200

