

## 兔抗 CHSY1 多克隆抗体

中文名称：兔抗 CHSY1 多克隆抗体

英文名称： Anti-CHSY1 rabbit polyclonal antibody

别名： CHSY; CSS1; TPBS; ChSy-1

抗原： CHSY1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	his gene encodes a member of the chondroitin N-acetylglucosaminyltransferase family. These enzymes possess dual glucuronyltransferase and galactosaminyltransferase activity and play critical roles in the biosynthesis of chondroitin sulfate, a glycosaminoglycan involved in many biological processes including cell proliferation and morphogenesis. Decreased expression of this gene may play a role in colorectal cancer, and mutations in this gene are a cause of temtamy preaxial brachydactyly syndrome.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CHSY1

<b>Immunogen:</b>	Synthesized peptide derived from internal of human CHSY 1.
<b>Full name:</b>	chondroitin sulfate synthase 1
<b>Synonyms :</b>	CHSY; CSS1; TPBS; ChSy-1
<b>SwissProt:</b>	Q86X52
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	92 kDa
<b>WB Positive control:</b>	HepG2 cells and K562 cells lysates
<b>WB Recommended dilution:</b>	500-3000



