

## 兔抗 AMOTL1 多克隆抗体

中文名称：兔抗 AMOTL1 多克隆抗体

英文名称：Anti-AMOTL1 rabbit polyclonal antibody

级别：JEAP

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a peripheral membrane protein that is a component of tight junctions or TJs. TJs form an apical junctional structure and act to control paracellular permeability and maintain cell polarity. This protein is related to angiotensin, an angiotensin binding protein that regulates endothelial cell migration and capillary formation. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	AMOTL1
<b>Immunogen:</b>	Synthetic peptide of human AMOTL1
<b>Full name:</b>	angiotensin like 1

<b>Synonyms:</b>	JEAP
<b>SwissProt:</b>	Q8IY63
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	107 kDa
<b>WB Positive control:</b>	A172 cell lysate
<b>WB Recommended dilution:</b>	200-1000

