

## 兔抗 CHUK(Phospho-Ser180/181) 多克隆抗体

中文名称：兔抗 CHUK(Phospho-Ser180/181) 多克隆抗体

英文名称： Anti-CHUK(Phospho-Ser180/181) rabbit polyclonal antibody

别名： IKK1; IKKA; IKBKA; TCF16; NFKBIKA; IKK-alpha

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： CHUK(Phospho-Ser180/181)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the serine/threonine protein kinase family. The encoded protein, a component of a cytokine-activated protein complex that is an inhibitor of the essential transcription factor NF-kappa-B complex, phosphorylates sites that trigger the degradation of the inhibitor via the ubiquitination pathway, thereby activating the transcription factor.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CHUK(Phospho-Ser180/181)

<b>Immunogen:</b>	Peptide sequence around phosphorylation site of serine 180/181 (C-T-S(p)-F-V) derived from Human IKK- alpha/ beta.
<b>Full name:</b>	conserved helix-loop-helix ubiquitous kinase
<b>Synonyms :</b>	IKK1; IKKA; IKBKA; TCF16; NFKB1A; IKK-alpha
<b>SwissProt:</b>	O15111
<b>WB Predicted band size:</b>	85 kDa
<b>WB Positive control:</b>	HepG2 cells lysates treated with TNF
<b>WB Recommended dilution:</b>	500-1000

