

兔抗 AMPD2 多克隆抗体

中文名称：兔抗 AMPD2 多克隆抗体

英文名称： Anti-AMPD2 rabbit polyclonal antibody

别名： PCH9; SPG63

抗原： AMPD2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is important in purine metabolism by converting AMP to IMP. The encoded protein, which acts as a homotetramer, is one of three AMP deaminases found in mammals. Several transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	AMPD2
Immunogen:	Synthesized peptide derived from internal of human AMPD2.
Full name:	adenosine monophosphate deaminase 2

Synonyms :	PCH9; SPG63
SwissProt:	Q01433
WB Predicted band size:	101 kDa
WB Positive control:	Hela cells, NIH/3T3 cells and COLO205 cells lysates
WB Recommended dilution:	500-3000

