

兔抗 AMPH 多克隆抗体

中文名称：兔抗 AMPH 多克隆抗体

英文名称： Anti-AMPH rabbit polyclonal antibody

别名： AMPH1

抗原： AMPH

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein associated with the cytoplasmic surface of synaptic vesicles. A subset of patients with stiff-man syndrome who were also affected by breast cancer are positive for autoantibodies against this protein. Alternate splicing of this gene results in two transcript variants encoding different isoforms. Additional splice variants have been described, but their full length sequences have not been determined. A pseudogene of this gene is found on chromosome 11.
Applications:	ELISA, WB, IHC
Name of antibody:	AMPH

Immunogen:	Fusion protein of human AMPH
Full name:	amphiphysin
Synonyms :	AMPH1
SwissProt:	P49418
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human ovarian cancer and Human lung cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	76 kDa
WB Positive control:	Mouse lung tissue
WB Recommended dilution:	500-2000



