

## 兔抗 AMY2A 多克隆抗体

中文名称：兔抗 AMY2A 多克隆抗体

英文名称： Anti-AMY2A rabbit polyclonal antibody

别名： PA; AMY2; AMY2B

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： AMY2A

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Amylases are secreted proteins that hydrolyze 1,4-alpha-glucoside bonds in oligosaccharides and polysaccharides, and thus catalyze the first step in digestion of dietary starch and glycogen. The human genome has a cluster of several amylase genes that are expressed at high levels in either salivary gland or pancreas. This gene encodes an amylase isoenzyme produced by the pancreas.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	AMY2A
<b>Immunogen:</b>	Fusion protein of human AMY2A
<b>Full name:</b>	amylase, alpha 2A (pancreatic)

<b>Synonyms :</b>	PA; AMY2; AMY2B
<b>SwissProt:</b>	P04746
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human lung cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	58 kDa
<b>WB Positive control:</b>	Mouse pancreas tissue
<b>WB Recommended dilution:</b>	1000-5000



