

## 兔抗 ANAPC1 多克隆抗体

中文名称：兔抗 ANAPC1 多克隆抗体

英文名称： Anti-ANAPC1 rabbit polyclonal antibody

别名： APC1; MCPR; TSG24

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ANAPC1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a subunit of the anaphase-promoting complex. This complex is an E3 ubiquitin ligase that regulates progression through the metaphase to anaphase portion of the cell cycle by ubiquitinating proteins which targets them for degradation. Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ANAPC1

<b>Immunogen:</b>	Synthetic peptide of human ANAPC1
<b>Full name:</b>	anaphase promoting complex subunit 1
<b>Synonyms:</b>	APC1; MCPR; TSG24
<b>SwissProt:</b>	Q9H1A4
<b>IHC positive control:</b>	Human gastric cancer and human liver cancer
<b>IHC Recommend dilution:</b>	15-50

