

兔抗 CIB2 多克隆抗体

中文名称：兔抗 CIB2 多克隆抗体

英文名称： Anti-CIB2 rabbit polyclonal antibody

别名： KIP2; USH1J; DFNB48

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： CIB2

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is similar to that of KIP/ CIB, calcineurin B, and calmodulin. The encoded protein is a calcium-binding regulatory protein that interacts with DN A-dependent protein kinase catalytic subunits (DNA-PKcs), and it is involved in photoreceptor cell maintenance. Mutations in this gene cause deafness, autosomal recessive, 48 (DFNB48), and also Usher syndrome 1J (USH1J).
Applications:	IHC
Name of antibody:	CIB2
Immunogen:	Synthesized peptide derived from internal of human CIB2.
Full name:	calcium and integrin binding family member 2

Synonyms :	KIP2; USH1J; DFNB48
SwissProt:	O75838
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

