

## 兔抗 CILK1 多克隆抗体

中文名称：兔抗 CILK1 多克隆抗体

英文名称： Anti-CILK1 rabbit polyclonal antibody

别名： ciliogenesis associated kinase 1; ECO; ICK; MRK; L

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CILK1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Eukaryotic protein kinases are enzymes that belong to a very extensive family of proteins which share a conserved catalytic core common with both serine/threonine and tyrosine protein kinases. This gene encodes an intestinal serine/threonine kinase harboring a dual phosphorylation site found in mitogen-activating protein (MAP) kinases. The protein localizes to the intestinal crypt region and is thought to be important in intestinal epithelial cell proliferation and differentiation. Alternative splicing has been observed at this locus and two variants, encoding the same isoform, have been identified. [provided by RefSeq, Jul 2008]

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CILK1
<b>Immunogen:</b>	Fusion protein of human CILK1
<b>Full name:</b>	ciliogenesis associated kinase 1
<b>Synonyms:</b>	ECO; ICK; MRK; LCK2; EJM10
<b>SwissProt:</b>	Q9UPZ9
<b>IHC positive control:</b>	Human thyroid cancer and Human brain
<b>IHC Recommend dilution:</b>	25-100



