

兔抗 ANGPTL2 多克隆抗体

中文名称：兔抗 ANGPTL2 多克隆抗体

英文名称：Anti-ANGPTL2 rabbit polyclonal antibody

别名：angiopoietin like 2; ARP2; HARP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ANGPTL2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Angiopoietins are members of the vascular endothelial growth factor family and the only known growth factors largely specific for vascular endothelium. Angiopoietin-1, angiopoietin-2, and angiopoietin-4 participate in the formation of blood vessels. A NGPTL2 protein is a secreted glycoprotein with homology to the angiopoietins and may exert a function on endothelial cells through autocrine or paracrine action.
Applications:	ELISA, IHC
Name of antibody:	ANGPTL2

Immunogen:	Synthetic peptide of human ANGPTL2
Full name:	angiopoietin like 2
Synonyms:	ARP2; HARP
SwissProt:	Q9UKU9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and Human esophagus cancer
IHC Recommend dilution:	50-100



