

## 兔抗 ANKRD10 多克隆抗体

中文名称：兔抗 ANKRD10 多克隆抗体

英文名称：Anti-ANKRD10 rabbit polyclonal antibody

别名：ankyrin repeat domain 10

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ANKRD10

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

ANKRD10 (Ankyrin repeat domain-containing protein 10) is a 420 amino acid protein that contains four Ankyrin repeats and is expressed as two alternatively spliced isoforms. The gene encoding ANKRD10 maps to chromosome 13. Comprising nearly 4% of the human genome, chromosome 13 contains around 114 million base pairs and encodes over 400 genes. Chromosome 13 houses key tumor suppressor genes, including BRCA2 and RB1, which are associated with breast cancer susceptibility and retinoblastoma, respectively. Trisomy 13, also known as Patau syndrome, is

	deadly and the few who survive past one year suffer from permanent neurologic defects, difficulty eating and vulnerability to serious respiratory infections.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ANKRD10
<b>Immunogen:</b>	Fusion protein of human ANKRD10
<b>Full name:</b>	ankyrin repeat domain 10
<b>SwissProt:</b>	Q9NXR5
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	200-300

