

兔抗 ANKRD54 多克隆抗体

中文名称：兔抗 ANKRD54 多克隆抗体

英文名称： Anti-ANKRD54 rabbit polyclonal antibody

别名： LIAR

相关类别： 一抗

抗原： ANKRD54

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Plays an important role in regulating intracellular signaling events associated with erythroid terminal differentiation. Interacts (via ankyrin repeat region) with LYN (via SH3-domain) in an activation-independent status of LYN By similarity. Forms a multiprotein complex with LYN and HCLS1 By similarity. Interacts with TSN2, VAV1, DBNL AND LASP1.
Applications:	ELISA, IHC
Name of antibody:	ANKRD54
Immunogen:	Fusion protein of human ANKRD54
Full name:	ankyrin repeat domain 54

Synonyms :	LIAR
SwissProt:	Q6NXT1
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human thyroid cancer
IHC Recommend dilution:	25-100

