

兔抗 ANO9 多克隆抗体

中文名称：兔抗 ANO9 多克隆抗体

英文名称： Anti-ANO9 rabbit polyclonal antibody

别名： PIG5; TP53I5; TMEM16J

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： ANO9

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Does not exhibit calcium-activated chloride channel (CaCC) activity. Can inhibit the activity of ANO1.
Applications:	IHC, IF
Name of antibody:	ANO9
Immunogen:	Synthesized peptide derived from internal of human TM16J.
Full name:	anoctamin 9
Synonyms :	PIG5; TP53I5; TMEM16J
SwissProt:	A1A5B4
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
IF positive control:	HepG2 cells

IF Recommend dilution:	100-500
------------------------	---------



