

兔抗 CLIC3 多克隆抗体

中文名称：兔抗 CLIC3 多克隆抗体

英文名称： Anti-CLIC3 rabbit polyclonal antibody

储 存： 冷冻（-20℃） 避光

抗 原： CLIC3

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Chloride channels are a diverse group of proteins that regulate fundamental cellular processes including stabilization of cell membrane potential, transepithelial transport, maintenance of intracellular pH, and regulation of cell volume. Chloride intracellular channel 3 is a member of the p64 family and is predominantly localized in the nucleus and stimulates chloride ion channel activity. In addition, this protein may participate in cellular growth control, based on its association with ERK7, a member of the MAP kinase family.
Applications:	WB
Name of antibody:	CLIC3
Immunogen:	Synthesized peptide derived from internal of human CLIC3.
Full name:	chloride intracellular channel 3

Synonyms :	0
SwissProt:	O95833
WB Predicted band size:	27 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

