

兔抗 APLP2(Ab-755) 多克隆抗体

中文名称：兔抗 APLP2(Ab-755) 多克隆抗体

英文名称：Anti-APLP2(Ab-755) rabbit polyclonal antibody

别名：APPH; APPL2; CDEBP; APLP-2

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：APLP2(Ab-755)

宿主：Rabbit

反应种属：Human Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in the regulation of hemostasis. The soluble form may have inhibitory properties towards coagulation factors. May interact with cellular G-prote in signaling pathways. May bind to the DNA 5'-GTCA CATG-3'(CDEI box). Inhibits trypsin, chymotrypsin, plas min, factor XIA and plasma and glandular kallikrein. Modulates the Cu/Zn nitric oxide-catalyzed autodegra dation of GPC1 heparan sulfate side chains in fibrobl asts By similarity.
Applications:	WB
Name of antibody:	APLP2(Ab-755)

Immunogen:	Synthesized non-phosphopeptide derived from human APLP2 around the phosphorylation site of tyrosine 755 (P-T-Y(p)-K-Y).
Full name:	amyloid beta (A4) precursor-like protein 2
Synonyms :	APPH; APPL2; CDEBP; APLP-2
SwissProt:	Q06481
WB Predicted band size:	87 kDa
WB Positive control:	Rat brain tissue lysate
WB Recommended dilution:	500-3000

