

兔抗 APOA4 多克隆抗体

中文名称：兔抗 APOA4 多克隆抗体

英文名称：Anti-APOA4 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

抗 原：APOA4

宿 主：Rabbit

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Apolipoprotein (apo) A-IV gene contains 3 exons separated by two introns. A sequence polymorphism has been identified in the 3'UTR of the third exon. The primary translation product is a 396-residue preprotein which after proteolytic processing is secreted its primary site of synthesis, the intestine, in association with chylomicron particles. Although its precise function is not known, apo A-IV is a potent activator of lecithin-cholesterol acyl transferase in vitro.
Applications:	ELISA, WB, IHC
Name of antibody:	APOA4
Immunogen:	Fusion protein of human APOA4
Full name:	apolipoprotein A-IV

SwissProt:	P06727
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and human cervical cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	45 kDa
WB Positive control:	Human serum
WB Recommended dilution:	200-1000



