

## 兔抗 CLIP4 多克隆抗体

中文名称：兔抗 CLIP4 多克隆抗体

英文名称： Anti-CLIP4 rabbit polyclonal antibody

别名： RSNL2

抗原： CLIP4

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	4 isoforms of the human protein are produced by alternative splicing.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CLIP4
<b>Immunogen:</b>	Synthesized peptide derived from internal of human CLIP4.
<b>Full name:</b>	CAP-GLY domain containing linker protein family, member 4
<b>Synonyms :</b>	RSNL2
<b>SwissProt:</b>	Q8N3C7
<b>WB Predicted band size:</b>	76 kDa
<b>WB Positive control:</b>	HepG2 cells and Jurkat cells lysates
<b>WB Recommended dilution:</b>	500-3000

