

兔抗 APOE 多克隆抗体

- 中文名称：兔抗 APOE 多克隆抗体
- 英文名称：Anti-APOE rabbit polyclonal antibody
- 别名：AD2; LPG; LDLCQ5
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：APOE
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：Unconjugate

技术规格

Background:

Chylomicron remnants and very low density lipoprotein (VLDL) remnants are rapidly removed from the circulation by receptor-mediated endocytosis in the liver. Apolipoprotein E, a main apoprotein of the chylomicron, binds to a specific receptor on liver cells and peripheral cells. ApoE is essential for the normal catabolism of triglyceride-rich lipoprotein constituents. The APOE gene is mapped to chromosome 19 in a cluster with APOC1 and APOC2. Defects in apolipoprotein E result in familial dysbetalipoproteinemia, or type III hyperlipoproteinemia (HLP III), in which increased plasma cholesterol and trig

	lycerides are the consequence of impaired clearance of chylomicron and VLDL remnants
Applications:	WB
Name of antibody:	APOE
Immunogen:	Synthesized peptide derived from C-terminal of human APOE.
Full name:	apolipoprotein E
Synonyms :	AD2; LPG; LDLCQ5
SwissProt:	P02649
WB Predicted band size:	36 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

