

## 兔抗 APOL2 多克隆抗体

中文名称： 兔抗 APOL2 多克隆抗体

英文名称： Anti-APOL2 rabbit polyclonal antibody

别名： APOL3; APOL-II

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： APOL2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. Two transcript variants encoding the same protein have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	APOL2
<b>Immunogen:</b>	Synthetic peptide of human APOL2
<b>Full name:</b>	apolipoprotein L, 2
<b>Synonyms :</b>	APOL3; APOL-II

<b>SwissProt:</b>	Q9BQE5
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	37 kDa
<b>WB Positive control:</b>	Human fetal liver tissue and huvec cells
<b>WB Recommended dilution:</b>	500-2000



